Injection and Blood Sampling Practice Model Type I

with storage case

Outline

A series of techniques such as confirmation of injection site, insertion hypodermic needle and injection of medicine are possible with this model. The pressurized blood vessel can be palpated, and needle insertion provides a feeling similar to that offered by a real human arm. Blood vessel tube is durable enough to be used in numerous practice injections.



Needle injection has same feel as with a human arm.



Replaceable outer skin and blood vessel tube.




Features

- Made of a special silicone rubber, this lifelike model limb features synthetic skin that is remarkably similar to human skin.
- The sensation of injecting a hypodermic needle into this model's simulated skin and blood vessels is identical to that encountered when performing an injection on an actual patient. Practice of both medicine injection and blood sampling using this model provides trainees with a remarkably realistic experience.
- Resilient synthetic skin conceals injection marks well, even after numerous practice injections. Injection points are highly resistant to tears and other damage.
- The artificial blood features the same viscosity and color as those of actual human blood. Since it is entirely water soluble, spills are cleaned easily without the worry of staining.
- The model's lifelike circulation operates through a simple manual process. This rugged, serviceable model incorporates extremely reliable mechanisms for trouble-free service over the long term.
- This simply designed model enables easy exchange of both the synthetic skin and veins.

Components

- ① Arm skin
- Arm musculature
- ③ Blood vessel tube
- ④ Stand
- ⑤ Irrigator
- 6 Simulated blood
- ⑦ Elbow rest



Accessory Set LM-028A

1pc of arm skin, 5pcs of blood vessel tube, 500ml of simulated blood

Sold Separately



Specifications

·			
Arm	Approx.59(L) $ imes$ 9(W) $ imes$ 9(H)cm	Approx.2kg	
Stand	Approx.60(L) \times 16(W) \times 67(H)cm	Approx.1.5kg	
Storage case	Approx.67(L) \times 30(W) \times 21(H)cm	Approx.4kg	

Components

Componente			
Arm skin	1		
Arm musculature	1		
Blood vessel tube	6		
Base and stand	1		
Irrigator (500ml)	2		
Simulated blood (500ml)	1		
Elbow rest	1		
Pinchcock	2		
Tube connector	2		
Baby powder	1		
Transfusion set	1		
Storage case	1		

Spares

opuloo				
LM-028B	Arm skin	1		
LM-028C	Blood vessel tube	10		
LM-028D	Simulated blood (500ml)	1		
LM028E	Tube connector	2		
LM-028S-T	Arm Skin (Thick type)	1		

38